




www.immunostep.com



# BASOSTEP

Determination of basophil activation upon allergen estimation by Flow Cytometry.




## PRINCIPLE

Using defined allergens, this test provides information concerning of IgE-mediated allergic disorders via the CD63 antigen as a marker of basophil activation.


## ADVANTAGES

1. Safer, since it does not expose the patient to any allergens. Allergy blood testing is the preferred test for infants and very young children.
2. More specific identification of basophils than other thest. CD203c marker included in the test is specific marker for basophils and mast cells.
3. Simple and fast, the kit contains positive and negative controls in liophilized tubes, for ease handling.
4. Versatil kit that offers 3 differents solutions (antibody combinations, prices) depending on the clinician needs.




# ASSAY PROCEDURE

## 1. SAMPLE PREPARATION




## 2. DEGRANULATION




We recommended testing several dilution of each Allergen (e.g. 1:100; 1:1000; 1:10.000; 1:100.000; 1:1.000.000)

## 3. LABELLING




## 4. LYSING AND FIXATION



## 5. ADQUISITION

Set on the cytometer to store only the events in the region RI as imagen right:



Acquire and store all RI events possible. It is recommended to acquire at a low or medium speed to avoid cell aggregates.

BSTP-100T	BSTP2-100T	BST3-100T
<b>STAINING REAGENT</b> <ul style="list-style-type: none"> <li>CD63 FITC</li> <li>CD203c PE</li> <li>HLA-DR PERCP</li> <li>CD123 APC</li> </ul>	<b>STAINING REAGENT</b> <ul style="list-style-type: none"> <li>CD63 FITC</li> <li>CD123 PE</li> <li>HLA-DR PERCP</li> </ul>	<b>STAINING REAGENT</b> <ul style="list-style-type: none"> <li>CD63 FITC</li> <li>CD123 PE</li> <li>HLA-DR APC</li> </ul>
<b>POSITIVE CONTROL</b> <ul style="list-style-type: none"> <li>fMLP</li> <li>Anti-Ig E</li> </ul>		
<b>STIMULATION BUFFER</b> Containing calcium, heparin and IL-3		
<b>Required apparatus:</b> Flow cytometer with 488 nm and 633 nm excitation wavelength	<b>Required apparatus:</b> Flow cytometer with 488 nm excitation wavelength	<b>Required apparatus:</b> Flow cytometer with 488 nm and 633 nm excitation wavelength

## REFERENCES:

Sainte-Laudy, J, et al. Analysis of membrane expression of the CD63 human basophil activation marker. Applications to allergologic diagnosis. Allerg. Immunol. Paris 26, 211-4 (1994)



Address: Avda. Universidad de Coimbra, s/n  
 Cancer Research Center (C.I.C.)  
 Campus Miguel de Unamuno  
 37007 Salamanca (Spain)  
 Tel. / Fax: (+34) 923 294 827  
 E-mail: info@immunostep.com  
 www.immunostep.com