

SARS-CoV-2 Multiplex Serological Assay

Controls provided

Compatible with Standard cytometers

No specific software needed

CE IVD

Evaluation of IgG+IgA+IgM antibodies response to 4 different SARS-CoV-2 viral proteins simultaneously.

RBD

Spike

Nucleocapsid

Mpro*

Exclusive combination of 4 viral antigens analyzed simultaneously



Wide Linearity Range



Near to 100% of sensitivity and specificity



Differentiating natural infection from vaccine response



Reliable and reproducible results



More results per sample (>12 plex)



*Patented Protein for immunological analysis. Under **CSIC** Patent licence. Assay for detection of Cysteine-like Protease (Mpro) of SARS-CoV-2" | EP 203824958.

IgG+IgA+IgM Multiplex Microsphere-based Assay

All the information, just one sample.

This bead-based assay by flow cytometry has demonstrated to provide a wide range of information analyzing just one sample in a short period: detecting 4 viral proteins and 3 different immunoglobulins, simultaneously.





Figure 1: Graphical representation of Anti-SARS-CoV-2 flow cytometry bead-based array procedure. Source: Internal Data.

Differentiating Immune Response by Virus and Vaccination



Facilitating the massive screening of Covid-19 patients to evaluate their immune response.

Pfizer-BioNTech COVID-19 vaccine (Comirnaty) and COVID-19 positive samples were assayed to evaluate IgG antibody response to RBD, S, N and Mpro viral proteins.

MFI levels vary from one sample to another, indicating the existence of different antibody response profiles depending on the viral protein analyzed.

Figure 2: SARS-CoV-2 multiplex assay vaccinated samples data compared with positive and negative samples. Source: Internal data.

Providing Outstanding Values of Sensitivity and Specificity

	Days	TP	FN	PPA	95% (C.I)	TN	FP	NPA	95% (C.I)	Total
IgG	<7	30	5	88%	75% 100%	293	5	98%	94% 100%	298
	7-15	47	3	94%	87% 100%					
	>15	25	1	96%	86% 100%					
IgA	<7	33	5	83%	75% 99%	296	12	96%	91% 100%	308
	7-15	46	10	80%	72% 93%					
	>15	28	1	97%	87% 100%					
IgM	<7	29	7	81%	64% 97%	271	21	92%	88% 97%	292
	7-15	47	9	90%	66% 100%					
	>15	20	5	80%	57% 100%					

Table 1. Diagnostic sensitivity and specificity. Plasma specimens collected from patients with confirmed COVID-19 PCR positive results were tested with SARS-CoV-2 Multiplex IgG+IgA+IgM Assay. Negative percent agreement (NPA) was determined by using specimens collected prior to December 2019. TP, true positive; FN, false negative, PPA, positive percent agreement; TN, true negative; FP, false positive.



+34 923 29 48 27 | info@immunostep.com | www.sars-cov-2-test.immunostep.com



Proyecto subvencionado por el Ministerio de Industria, Comercio y Turismo, como parte del programa de ayudas para la fabricación de material sanitario por la crisis de la COVID-19 (2020).