


The diagram illustrates a 10-step protocol for the detection of IgG antibody in a sample. The process begins with the addition of 50ul of capture beads to a sample tube containing 50ul of sample. This is followed by a 30-minute room temperature (RT) orbital shake. Then, 1ml of wash buffer is added, and the mixture is centrifuged at 2500xg for 5 minutes. The supernatant is removed, and the beads are resuspended in 200ul of wash buffer. This is followed by another 30-minute RT orbital shake. The beads are then washed again by centrifuging at 2500xg for 5 minutes and removing the supernatant. Finally, 10ul of detector antibody is added, and the mixture is centrifuged at 2500xg for 5 minutes. The beads are resuspended in 200ul of wash buffer, and the mixture is centrifuged at 2500xg for 5 minutes. The final step is the detection of the antibody using a colorimetric assay.

ANNEX III. Serial dilution (1:2) of the standard for the construction of the standard curve



ANEXO IV. Interpretation of results

