

REF

MITO-100T



100 test



RUO

1. PRODUCT DESCRIPTION

Membrane potential ($\Delta\Psi$) is generated and maintained by concentration gradients of ions such as sodium, potassium, chloride, and hydrogen.

Mitochondrial $\Delta\Psi$ drives the accumulation in mitochondria of cationic dyes such as cyanines, and the mitochondrial $\Delta\Psi$ is reduced when energy metabolism is disrupted, notably in apoptosis. Changes in the mitochondrial $\Delta\Psi$ have been described during necrosis, cell cycle and apoptosis.

Mitochondrial uptake of dye is a possible source of fluorescence variance. Flow cytometry can be used to estimate membrane potential in eukaryotic cells. Methods using cyanines dyes can detect changes in $\Delta\Psi$.



Immunostep MitoStep uses a cationic dye DiICl(5) (1,1',3,3,3'-hexamethylindodicarbo-cyanine iodide) for the study of mitochondrial $\Delta\Psi$. During the apoptosis occurs depolarization of the membrane and as a result there is an increase in cells with less DiICl(5) fluorescence.

MitoStep has been optimized for use in flow cytometry, cells stained with DiICl(5) are excite using air-cooled Helium-Neon laser emitting at 633nm, cells DiICl(5) positives emitted at 658 nm, DiICl(5) mean intensity of fluorescence decreases when cells are treated with reagents that induce apoptosis or reagents that disrupt $\Delta\Psi$ mitochondrial.

Storage buffer: DiICl(5), 500 μ l of 10 μ M in DMSO.
Storage conditions: Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact our Technical Services. tech@immunostep.com

2. RECOMMENDATIONS AND WARNINGS

DMSO is a potentially toxic. It is recommended that the user wear protective clothing, gloves, and eye/face protection in order to avoid contact with the skin and eyes.

3. PROTOCOL

Staining cells protocol with DiICl(5)

1. Harvest the cells after the apoptosis induction or treatment with a disrupt membrane potential reagent and wash in temperate phosphate-buffered saline (PBS).
 2. Wash cells twice with temperate PBS and resuspend cells in temperate phosphate-buffered saline (PBS) at a concentration 1 x 10⁶ cells/ml.
 3. Add 5 μ l of 10 μ M DiICl(5).
 4. Incubate the cells at 37 °C, 5% CO₂, for 15 minutes.
 5. After incubation period, add 400 μ l of PBS to each tube. Analyze by flow cytometry.
- Please, refer to www.immunostep.com technical support for more information.




Figure 1. Jurkat cells (T-cell leukemia, human) treated with 6 μ M camptothecin for four hours (bottom panel) or untreated (top panel).

4. WARRANTY






Warranted only to conform to the quantity and contents stated on the label or in the product labelling at the time of delivery to the customer.

Immunostep disclaims hereby other warranties. Immunostep's sole liability is limited to either the replacement of the products or refund of the purchase price.

5. REFERENCES

1. Howard M. Shapiro. Membrane Potential Estimation by Flow Cytometry. Methods 279 (2000).
21, 271-

6. EXPLANATION OF SYMBOLS

	Form
<div>REF</div>	Catalog reference
	Contains sufficient for <n> test
	Quantity per test
	Regulatory Status
<div>RUO</div>	Research Use Only
	Manufacturer

7. MANUFACTURED BY: IMMUNOSTEP S.L.



Address: Avda. Universidad de Coimbra, s/n
Cancer Research Center (C.I.C)
Campus de Unamuno
37007 Salamanca (Spain)
Telf./fax: (+34) 923 294 827
E-mail: info@immunostep.com
www.immunostep.com