

Anti-Human CD4/CD8 (HP2-6/143-44)

Fluorochrome	Reference	Test
FITC/PE	4F18PEI-50T	50 tests



PRODUCT DESCRIPTION

Clone: HP2-6, 143-44

Isotype: Mouse IgG2a/ Mouse IgG1

Tested application: flow cytometry

Immunogen: The anti-CD4 monoclonal antibody derives from T cells from leukemic HPB-ALL. The anti-CD8 monoclonal antibody derives from human T cells.

Species reactivity: Human

Storage instruction: store in the dark at 2-8 °C

Storage buffer: aqueous buffered solution containing protein stabilizer and 0.09% sodium azide (NaN₃).

Recommended usage: Immunostep's CD4/CD8, is a two-color direct immunofluorescence reagent for enumerating percentages of mature human helper/inducer (CD4+) and suppressor/cytotoxic (CD8+) lymphocytes in erythrocyte-lysed whole blood (LWB). The helper/suppressor ratio (CD4+/CD8+) may also be determined.

Presentation: liquid

Source: Supernatant proceeding from an in vitro cell culture of a cell hybridoma.

Purification: Affinity chromatography

ANTIGEN DETAILS

Large description: Percentages of CD4+ and CD8+ lymphocytes are used in monitoring the immune status of patients with immune deficiency diseases, autoimmune diseases, or immune reactions.

The relative percentage of the CD4 + subset is depressed and the relative percentage of the CD8 +subset is elevated in many patients with congenital or acquired immune deficiencies 1 such as severe combined immunodeficiency (SCID) and acquired immunodeficiency syndrome (AIDS).

The percentage of suppressor/cytotoxic lymphocytes can be outside the normal reference range in some autoimmune diseases 4 and in certain immune reactions such as acute graft-versus-host disease (GVHD) and transplant rejection. The relative percentage of the CD8 + lymphocyte population may often be decreased in active systemic lupus erythematosus (SLE) but can also be increased in SLE patients undergoing steroid therapy.

The CD4 + /CD8 + (helper/suppressor) lymphocyte ratio, quantified as the ratio of CD4 fluorescein isothiocyanate (FITC)-positive lymphocytes to CD8 phycoerythrin (PE)-positive lymphocytes, has been used to evaluate the immune status of patients with, or suspected of developing, autoimmune disorders or immune deficiencies. In many cases, the relative percentages of helper lymphocytes decline and suppressor lymphocytes increase in immune deficiency states. These states may also be marked by T-cell lymphopenia. In addition, the ratio has been used to monitor bone marrowtransplant patients for onset of acute GVHD. While a useful indicator, the CD4 + /CD8 + (helper/ suppressor) lymphocyte ratio has specific limitations. (1-2)

Please, refer to www.immunostep.com technical support for more information.

INSTRUCTIONS FOR USE

Store Diluting RBCX10 Lysis Solution:

1. Dilute the 10X RBCX10 concentrate 1:10 with room temperature (20–25 °C) deionized water.
2. The prepared solution is stable for 1 month when stored in a glass or high-density polyethylene (HDPE) container at room temperature.

Performing Quality Control: In accordance with the College of American Pathologists (CAP) guidelines, we recommend running two levels of liquid control material (process control).

Controls should be run at least once each day that patient testing is performed. Use commercial controls providing established values for percent positive and absolute counts with each run to assess system performance. Immunostep recommends using CDChex Plus (Streck) normal and CD4 Low as process controls.

To perform quality control: Thoroughly mix the appropriate CD-Chex Plus (Streck) control, or equivalent process control. See the IFU for the control for detailed instructions.

1. Stain the control sample using Immunostep's CD4/CD8 antibody as described in the following section. The control sample should be processed like patient samples to monitor the ongoing performance of the entire analytic process.
2. Acquire the stained control sample on the flow cytometer.
3. Visually inspect the CD45 vs SSC dot plot. The lymphocyte population should appear as a bright, compact cluster with low SSC. Monocytes and granulocytes should also appear as distinct clusters. Do not proceed with analysis if populations are diffuse and there is little or no separation between clusters.
4. Verify that the results are within the values reported on the Assay Values sheet.

Staining the Cells:

Use care to protect the tubes from direct light. Perform the procedure at room temperature. See Precautions and Interfering Conditions.

1. For each patient sample, label a 12 × 75-mm tube with the sample identification number. For absolute counts, label an STEPCOUNT counting tube in place of the 12 × 75-mm tube. Note: Before use, verify that the STEPCOUNT bead pellet is intact and within the metal retainer at the bottom of the tube. If this is not the case, discard the STEPCOUNT counting tube and replace it with another.
2. Pipette 20 µL of Immunostep's CD4/CD8 antibody into the bottom of the tube. Note: If using an STEPCOUNT counting tube, pipette the reagent onto the side of the tube, just above the metal retainer, without touching the bead pellet.
3. Pipette 50 µL of well-mixed, anticoagulated whole blood into the bottom of the tube. Note: If using an STEPCOUNT counting tube, we recommend using the reverse pipetting technique to pipette the sample onto the side of the tube just above the metal retainer. Avoid smearing blood down the side of the tube. If whole blood remains on the side of the tube, it will not be stained with the reagent and can affect results.
4. Cap the tube and vortex gently to mix.
5. Incubate for 15–30 minutes in the dark at room temperature (20–25 °C).
6. Add 450 µL of IX RBCX10 lysis solution to the tube.
7. Cap the tube and vortex gently to mix.
8. Incubate for 15–30 minutes in the dark at room temperature (20–25 °C). The sample is now ready to

be analyzed on the flow cytometer. If samples will not be analyzed immediately after staining, store them in the dark at room temperature (20–25°C).

Acquiring the Samples:

1. Vortex the cells thoroughly at low speed. It is important to reduce aggregation before running samples on the flow cytometer. Note: If you are using a Loader, vortex tubes immediately before placing them into the Loader racks.
2. Install the tube on the cytometer and acquire the sample. Before acquiring samples, adjust the threshold to minimize debris and ensure that populations of interest are included.
3. Analyze the data using the appropriate cytometer-specific software. See the cytometer's IFU for more information.

Results: Results are reported as the percentage of positive cells per lymphocyte population or as the number of positive cells per microliter of blood (absolute count).

Calculating Absolute Counts: During analysis, the absolute number (cells/ μ L) of positive cells in the sample can be determined by comparing cellular events to bead events. If Immunostep software is used, absolute counts will be determined by the software. For manual data analysis, the absolute count of the cell population (A) can be calculated using the following equation:

$$A=X/Y \times N/V$$

Where:

- X is the number of positive cell events
- Y is the number of bead events
- N is the number of beads per test, which is found on the STEPCOUNT counting tubes foil pouch and can vary from lot to lot
- V is the sample volume (50 μ L)

Controls:

- Positive Controls: Use CD-Chex Plus (Streck) normal and CD4 Low to validate the staining protocol.
- Negative Controls: Include isotype controls and unstained cells to set proper gating strategies.

Material Not Supplied:

- STEPCOUNT counting tubes (ref. 1399991218)
- RBCX10 lysis solution (ref. RBCX10-50ML)
- CD-Chex Plus (Streck) normal and CD4 Low

REFERENCES

1. Giorgi JV. Characterization of T lymphocyte subset alterations by flow cytometry in HIV disease. *Ann N Y Acad Sci.* 1993 Mar 20;677:126-37. doi: 10.1111/j.1749-6632.1993.tb38771.x. PMID: 8494202.
2. Gratama JW, Naipal A, Olijans P, Zwaan FE, Verdonck LF, de Witte T, Vossen JM, Bolhuis RL, de Gast GC, Jansen J. T lymphocyte repopulation and differentiation after bone marrow transplantation. Early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. *Blood.* 1984 Jun;63(6):1416-23. PMID: 6372896.

MANUFACTURER BY:



Immunostep S.L
Avda. Universidad de Coimbra, s/n Cancer Research Center (CIC) Campus Miguel de Unamuno 37007 Salamanca (Spain)
Tel. (+34) 923 294 827
www.immunostep.com

SURVEILLANCE AND NOTIFICATION

In accordance with Annex I, Section 20.4.1(n) of Regulation (EU) 2017/746, the user is obliged to report any serious incident related to the use of the product.

- **To the Manufacturer:** Please contact our Surveillance Department at vigilancia@immunostep.com
- **To the Competent Authority:** Report via the official channels of the Member State.

EXPLANATION OF SYMBOLS



CE labeling



In vitro diagnostic device



Manufacturer



Pay attention to