

Wash Buffer (10X)

1. INTRODUCTION



PBS 10% BSA.

Reference: IMS0515

Volume: 20 mL

Storage: 2-8°C

pH-value: 7.4 ± 0.2

Preservative: contains 10 % albumin in 10mM sodium phosphate, 150 mM NaCl, pH 7.4, contains KATHON anti-microbial agent.

Expiry date when stored unopened: please refer to the label on the bottle.

FOR RESEARCH USE ONLY, NOT FOR DIAGNOSTIC USE

2. INSTRUCTIONS FOR USE

- Before use, temper the buffer between +18°C y +24°C (room temperature) for 30 minutes.
- If crystallization is observed in the concentrated buffer during storage, warm to 37°C and properly shake before diluting.
- Dilute the contents of the 10X assay buffer to 1X (PBS) in PBS 1X, for use flow cytometry assays (Exosomes and SARS-CoV-2 Immunostep Assays). To carry out the dilution, the amount required for the assays is removed from the concentrate bottle and diluted 1:10 in PBS 1X.

3. RELATED PRODUCTS

- Covid Multiplex (Ref. IMS0509)
- Exosomes (Kits and Capture Beads)
- HemoStep (Ref. IMS0212)
- Fusion Proteins

4. MANUFACTURED BY

IMMUNOSTEP S.L.

Address: Avda. Universidad de Coimbra, s/n
Cancer Research Center (C.I.C)
Campus de Unamuno
37007 Salamanca (Spain)
Telf./fax: (+34) 923 294 827
E-mail: info@immunostep.com
www.immunostep.com